



# DENTASEPT®

## Aspiration AF+

### INDICATIONS

Cleaning and disinfection of suction systems in a dental environment. Compatible with different suction systems. Does not produce foam.

### ADVANTAGES

- ✓ Non-foaming
- ✓ Triple Action - cleaning, descaling and disinfecting
- ✓ Proven efficiency on biofilm according to ISO/TS 15883 – 5 : 2006
- ✓ Very economical 2% dilution: 5L of concentrated solution = 250L of ready-to-use solution
- ✓ Effective within 15 minutes - active against all enveloped viruses<sup>1</sup>
- ✓ Not classified as dangerous for the environment when diluted at 2%
- ✓ Class IIa medical device (Directive 93/42/EEC)

<sup>1</sup> The Robert Koch Institute (RKI) recommends the use of disinfectants with a limited virucidal activity to stop the infection chain. According to RKI, a limited virucidal activity (effective on enveloped viruses) means an action on Vaccinia virus and BVDV.

### MICROBIOLOGICAL PROPERTIES AT 20°C

Bactericide: EN 1040, EN 13727, EN 14561 – 15 minutes  
 Yeasticide: EN 1275, EN 13624, EN 14562 – 15 minutes  
 Tuberculocide: EN 14348, EN 14563 – 60 minutes  
 Active according to EN 14476 on Herpesvirus, BVDV (Surrogate of HCV), Vaccinia virus – 5 minutes.

### COMPOSITION

N-(3-aminopropyl)-N-dodecylpropane-1,3-diamine, association of surfactants, sequestering agent, colorant and excipients.

### PACKAGING & ACCESSORIES

1 litre dosing bottle  
 5 litre container with a dosing pump  
 Dentasept® ASP System (424050).



(424050)



ASPIRATION

# DENTASEPT® Aspiration AF+

## PROTOCOL FOR USE

Cleaning and disinfection of suction systems in a dental environment. Compatible with different suction systems. Does not produce foam.

### BEFORE USE

Wear suitable gloves.

### HOW TO USE

At the end of the day, prepare 2 litres of diluted solution with the Dentasept®ASP System.

- 1 Dilute at 2% (40 ml in 2 litres of water).  
1 litre bottle: Fill the dispenser up to the 20 ml indicator, then repeat the procedure.  
5 litre container: 2 presses on the dosing pump.
- 2 Gently shake the Dentasept®ASP System to mix the solution.
- 3 Partially unscrew the black cap to let air into the Dentasept® ASP System. Attach the suction system tubes and aspirate the diluted solution.
- 4 Leave for a minimum of 15 minutes. The solution can be left in the suction system overnight.
- 5 To rinse the system, fill the Dentasept® ASP System with a minimum of 1 litre clean water and aspirate.

Any prolonged exposure to a high temperature (>35°C) may lead to a change in the colour without altering the antimicrobial performance.

### PRECAUTIONS

Dangerous: Strictly follow the precautions for use (developed in accordance with European regulations on the classification and labelling of chemicals). Before use, carefully read the label and product information. For professional use only. Class IIa medical device (in accordance with amended Directive 93/42/EEC on medical devices). Store between +5 and +35°C.

