



MEDICAL DEVICES

DENTASEPT®

SH PRO

#### INDICATIONS

Foam detergent disinfectant for the treatment of medical devices in the dental environment.

#### ADVANTAGES

- ✓ New improved formulation
- ✓ Ready to use solution
- ✓ Formulated without alcohol, compatible with all synthetic materials, plastics etc
- ✓ Without perfume, colourant or allergy promoting substances
- ✓ Foamer spray limits aerosolization of the product
- ✓ Effective in just 1 minute against all enveloped viruses<sup>1</sup>
- ✓ Medical Device Class IIb (Directive 93/42/EEC modified)

<sup>1</sup> The Robert Koch Institute recommends the use of disinfectants claiming limited virucidal activity to interrupt the infection-chain. According to the RKI, limited virucidal activity (effective against enveloped viruses) corresponds to an activity against Vaccinia virus and BVDV.

#### MICROBIOLOGICAL PROPERTIES AT 20°C

Bactericide: EN 1040, EN 13727, EN 13697 - 2 minutes.

MRSA: EN 13727, EN 13697 - 2 minutes.

Yeasticide: EN 1275, EN 13624, EN 13697 - 5 minutes.

Active according to EN 14476 against Rotavirus, Herpes virus, PRV (surrogate of HBV), BVDV (surrogate of HCV), Vaccinia virus - 1 minute.

Active according to EN 14476 against polyomavirus SV40 - 5 minutes.

#### COMPOSITION

Didecyldimethylammonium chloride, excipients.

#### PACKAGING

750ml Foamer Spray.

5 litre Container.



MEDICAL DEVICES

# DENTASEPT® SH PRO

## PROTOCOL FOR USE

Foam detergent disinfectant for the treatment of medical devices in the dental environment.

### BEFORE USE

Wear appropriate gloves.

### HOW TO USE

- 1 Spray Dentasept® SH PRO on to a non-woven cloth or directly on to the area to be cleaned and disinfected.
- 2 Evenly distribute the product (+/- 6 ml/m<sup>2</sup> or approx. 10 sprays) with a non-woven cloth to ensure that the complete area is treated.
- 3 Observe the contact time specified to achieve the level of activity required.  
Rinsing is not necessary unless the medical device comes into contact with skin or mucous membranes.

### PRECAUTIONS

Dangerous - respect the precautions for use. Medical Device Class IIb, CE 0459. Before use, read the label and information concerning the product.

Store between +5°C and +35°C.

