



MEDICAL DEVICES



DENTASEPT®

SH^{PRO} Wipes

INDICATIONS

Ready for use wipes impregnated with a cleaning and disinfectant solution for the treatment of invasive medical devices in dental environments such as wiping reusable instruments (rotary hand pieces, turbines, contra-angles, etc.) during preliminary treatment and wrapping to prevent the drying out of the soiled areas whilst proceeding to the next steps, and the curing lamp bodies (non-immersible and non-invasive medical devices).

ADVANTAGES

- ✓ Formulated without alcohol, compatible with all synthetic materials, plastics etc
- ✓ Without perfume, colourant or allergy promoting substances
- ✓ Wipes made from 100% biodegradable viscose
- ✓ Size: 180 x 200mm
- ✓ Medical Device Class IIb (Directive 93/42/EEC as amended)

MICROBIOLOGICAL PROPERTIES AT 20°C

Bactericide and yeasticide wipes in 2 min, 20°C (EN 16615).
 Impregnation solution: Bactericide in 1 min, 20°C (EN 13727).
 Yeasticide in 2 min, 20°C (EN 13624).
 Active in accordance with EN 14476 in 1 min, 20°C on HIV-1, PRV (surrogate of HBV), BVDV (surrogate of HCV), Vaccinia virus, Rotavirus, Herpesvirus, RSV.
 Active in accordance with EN 14476 in 5 min, 20°C on Coronavirus, and in 15 min, 20°C, on Norovirus.

COMPOSITION

Impregnation solution: Didecyldimethylammonium chloride, excipients.

PACKAGING

Pack of 100 wipes.



MEDICAL DEVICES

DENTASEPT® SH® Wipes

PROTOCOL FOR USE

Wipes impregnated with a cleaning and disinfectant solution for the treatment of medical devices in the dental environment.

BEFORE USE

Wear appropriate gloves.

HOW TO USE

- 1 Remove wipe from pack.
- 2 Apply to the area of the medical device to be treated, using as many wipes as necessary.
- 3 Observe the contact time specified to achieve the level of activity required.*
- 4 Rinsing is not necessary unless the medical device comes into contact with skin or mucous membranes.

Close the packaging properly after each time it is opened. Do not reuse wipe. Disinfection can no longer be guaranteed. Scientific dossier (efficacy, compatibility, etc.) available on request.

PRECAUTIONS

Dangerous – respect the precautions for use (Drawn up according to the European rules in force regarding the classification and labelling of chemical products). Before every use, read the label and information concerning the product.

Class IIb medical device (Directive 93/42/EEC as amended). Reserved for professional use. Store between +5°C and +35°C.

