

LUXATORS LUXATOREN




3 mm	190/L1ST-ED	190/L1C-ED	190/3ST-ED	190/3C-ED
4 mm	190/L2ST-ED	190/L2C-ED	190/4ST-ED	190/4C-ED
4,5 mm			190/45ST-ED	190/45C-ED

The LEIBINGER ERGODESIGN


Luxators are used for severing the delicate fibres connecting the tooth root with the jaw bone (ligaments) and extending the tooth socket (alveolus) thus releasing the tooth from the bone and allowing it to be removed easily. This allows the tooth to be extracted gently, which is ideal preparation for implantation.

Mit den LEIBINGER ERGODESIGN

Luxatoren können die feinen Fasern, die die Zahnwurzel mit dem Kieferknochen verbinden (Ligamente) durchtrennen und das Zahnfach (Alveole) aufgeweitet werden, so dass der Zahn die Verbindung zum Knochen verliert und leicht entfernt werden kann. Somit erhalten Sie eine schonende Extraktion, ideal für die Vorbereitung einer Implantation.



NEW	2,5mm	190/25C	190/25ST
	3 mm	190/3C	190/3ST
	4 mm	190/4C	190/4ST
	5 mm	190/5C	190/5ST



3 mm	190/3CF	190/3STF
4 mm	190/4CF	190/4STF
5 mm	190/5CF	190/5STF

LUXATOREN LUXATORS