

KIDS CLEAR CROWN

C L E A R
P R I M A R Y
A N T E R I O R
C R O W N
F O R M S



OUR KIDS' CROWNS



Products qualified by doctors and patients worldwide



Extremely safe with Rings



Golden titanium is coated on ss-crown for esthetic choice



Metal-free, esthetic solution

Features



Aesthetic Excellence

Thin interproximal walls.
Automatically contour itself to match the natural shape.



Easy Handling

Sufficient strength.
Strip off easily leaving a smooth surface.



Highly Effective

High preservation rate.
Simplify composite work for pediatric anterior restorations.

Product Composition

Central & Lateral KIT

A1 Upper Right *6.0mm	A2 Upper Right *6.7mm	A1 Upper Left *6.0mm	A2 Upper Left *6.7mm
A3 Upper Right *7.4mm	A4 Upper Right *8.1mm	A3 Upper Left *7.4mm	A4 Upper Left *8.1mm
B1 Upper Right *4.3mm	B2 Upper Right *5.1mm	B1 Upper Left *4.3mm	B2 Upper Left *5.1mm
B3 Upper Right *5.9mm	B4 Upper Right *6.7mm	B3 Upper Left *5.9mm	B4 Upper Left *6.7mm

KIT DESCRIPTIONS

- **TYPES** | A - Central
B - Lateral
- **SIZE** | 1<2<3<4
* Mesial-Distal (mm)



Shinhung.Co.,Ltd.

450, Cheongpa-ro, Jung-gu, Seoul, Korea
Tel : 82-2-6366-2175 / Fax : 82-2-6366-2131
<http://www.shinhung.com>
E-mail : export@shinhung.co.kr

We reserve the rights to modify the design or specification without prior notice to allow for the improvement.