

KIDS CROWN FAMILY

Deciduous teeth

YES

NO

Aesthetic Needs

Affordable Quality

More Safety

For patients of all ages

KIDS CROWN

Our crowns are loved by doctors and patients all over the world!



KIDS GOLDEN CROWN

Titanium coated crowns now have a natural, lustrous gold colour



KIDS RING CROWN

Rings add extra safety to our existing crowns!



SS Crown for Permanent Teeth, appropriate for patients of all ages from small children to the elderly

PERMA CROWN



PERMA GOLDEN CROWN



Kids crowns are loved by doctors and patients all over the world!

Check out our crowns line up today!



KIDS CROWN

Our crowns are loved by doctors and patients all over the world!



KIDS RING CROWN

Rings add extra safety to our existing crowns!



KIDS GOLDEN CROWN

Titanium coated crowns now have a natural, lustrous gold colour



SHINHUNG

Shinhung, Co., Ltd.
450, Cheongpa-ro, Jung-gu, Seoul, Korea
Tel : 82-2-6366-2175 / Fax : 82-2-6366-2131
<http://www.kidscrown.com>
E-mail : export@shinhung.co.kr

We reserve the rights to modify the design or specification without prior notice to allow for the improvement.



Our Crowns are loved by doctors and patients all over the world!

KIDS CROWN



Kids Crown Anterior Kit/Refill

Kids Crown Posterior Kit/Refill

Features

- 1 Pre-Contoured Design**
Save valuable chair time with our optimized pre-contoured SSCs, perfectly designed for the Hall technique.
- 2 Transparent Ceramic Coating**
Reduce nickel allergy risk and achieve a natural look with our transparent, gloss-minimizing ceramic coating.
- 3 Strong and Flexible Occlusal Surface**
Soft enough to prevent perforation, yet durable for lasting performance.
- 4 Upper D Type with Shorter Mesiodistal Length**
Perfect fit ensured.

Composition

KIDS CROWN POSTERIOR KIT

- 96 EA (48 sizes X 2 EA per size)
- 48 EA (48 sizes X 1 EA per size)

KIDS CROWN ANTERIOR KIT

- 36 EA (36 sizes X 1 EA per size)

KIDS CROWN REFILL

- 5 EA (same size)



Our White Crowns provide children with a beautiful and confident smile.

KIDS ZIRCONIA CROWN



Central&Lateral/Anterior full Kit/Refill
Posterior Kit/Refill

Central&Lateral/Anterior full Try-in Kit/Refill
Posterior Try-in Kit/Refill

Features

- 1 High retention**
Sand-blasting maximizes cement-surface bond
- 2 Minimal Thickness**
Thinner by 0.2mm, reducing prep time
- 3 Natural color and translucency**
Similar to a natural tooth in color and translucency, which is ideal for primary teeth.
- 4 Wide Open Gingival Margin**
Designed for easier fitting and optimal cingulum preservation.
- 5 Narrow Type Design**
Narrow crown design is ideal for limited spaces.

Composition

Type			Zirconia(For Restoration)	TRY-IN(For Trial)
Anterior	Anterior Full	KIT	48 Crowns (48 sizes*1EA)	48 Crowns (48 sizes*1EA)
	Central & Lateral	KIT	24 Crowns (24 sizes*1EA)	24 Crowns (24 sizes*1EA)
		REFILL	2EA(Same size)	1EA
	Cuspid&Cuspid Narrow	REFILL	2EA(Same size)	1EA
Universal (Mandibular)	REFILL	2EA(Same size)	1EA	
Posterior	Regular Posterior	KIT	40 Crowns (40 sizes*1EA)	40 Crowns (40 sizes*1EA)
		REFILL	2EA(Same size)	1EA
	Narrow Posterior	KIT	40 Crowns (40 sizes*1EA)	40 Crowns (40 sizes*1EA)
		REFILL	2EA(Same size)	1EA