

# PERMA CROWN

STAINLESS STEEL PERMANENT MOLAR CROWN



## Perma Crown is Excellent for

- ✓ Permanent teeth in children or adults
- ✓ Restoration of immature molar enamel or dentin abnormalities (Adolescent patients)
- ✓ Temporary use (After implant fixture placement or Temporary prosthesis)
- ✓ Observation after endo treatment with a poor prognosis and possibility of implantation

# PERMA CROWN



## Kit Configurations

**Perma Crown KIT**

**24pcs** (1EA per size)

**Perma Crown Refill**

**Refill boxes** (5EA / 10EA)

## Components

\*\* mm = Mesial-Distal

